Rotary Evaporator

Highly efficient standard rotary evaporator with manual lift

RE202-A/212-A (basic) REV202M-A/212M-A (with vacuum controller)



Features

5~315 rpm rotation speed range Turning the encoder dial slowly increases or decreases the value by 1, turning it quickly changes the value by 10. Selectable rotation mode (forward, reverse, auto reverse) RPM display brightness can be adjusted in 8 levels.

Set inversion function

Glassware and bath can be set in either side, left or right, depending on user's dominant hand and installation location. Three units of Glassware B (with vertical condenser) can be installed in a standard fume hood.



Uniquely designed glass condenser

This prevents liquid stagnation and backflow improving durability of vacuum seal

When using ketone or ether solvents, standard vacuum seal swells. It is recommended to use PTFE vacuum seal.





Two types of baths

Option for water and oil bath.

Large capacity 5L bath with 240mm I.D. and a full complement of safety functions such as automatic overheat prevention and temp. upper limit difference.



Specifications

Model		WITHOUT VACUU	JM CONTROLLER	WITH VACUUM CONTROLLER			
		RE202-AWA/212-AWA (Glassware A)	RE202-BWA/212-BWA (Glassware B)	REV202M-AWA/212M-AWA (Glassware A)	REV202M-BWA/212M-BWA (Glassware B)		
	Operating ambient temp. range	5~35°C					
Performance *1	Speed range	5~315 rpm * ³					
renormance	Evaporation capacity	Up to 23 ml/min					
	Pressure setting range			0-101	3 hPa		
	RPM display		Digital disp	lay / Control knob			
Functions	Rotation mode		Forward / Rev	verse / Auto inversion			
	Spring-loaded jack	Manual	Manual balance (maximum height 200 mm, stepless regulation, one-touch lock)				
	Vacuum controller			VR102S, installed above jack h	andle with attachment bracket		
Configuration	Vacuum control solenoid valve			OVR10, installed in the	ne rear of stand base		
Comgulation	Rotary motor		DC brushle	ess (simple servo)			
	Condenser retention		Condenser bracket		Condenser bracket		
Safety functions	Rotary evaporator		otor overload protection, overvoltage, low voltage, rotation speed sensor error Short circuit in internal circuit, overcurrent protection, overvoltage protection				
	Vacuum controller			Communication error, Pressure sensor error, Memory err Leak error, High pressure error, Auto leak at error occurre			
	Cooling condenser	Double corrugated tube (cooling surface: 0.143 m ²)					
		Suction port: GL-14 (lower), Φ10 nozzle	Suction Port: GL-14 (upper), Φ10 nozzle	Suction port: GL-14 (lower), Φ10 nozzle	Suction Port: GL-14 (upper), Φ10 nozzle		
		Cooling port: GL14 (two places in lower part), two ϕ 10 nozzles					
	Compatible evaporation flask		50-2000ml. Use optional reducer to attach small flasks				
	Compatible receiving flask	100-2000 ml					
Standard	External dimensions *2 (W x D x H)	719 × 324 × 534	529 × 324 × 745	719 × 324 × 534	529 × 324 × 745		
	Overall dimensions *2 (including bath) (W x D x H)	744 × 365 × 534	554 × 365 × 745	744 × 365 × 534	554 × 365 × 745		
	Power rating 50/60Hz	RE202 and REV202M: 100-115V 1A <i>with plug</i> RE212 and REV212M: 200-230V single phase 1A <i>no plug, round terminal</i>					
	Weight	~10.0 kg ~10.5 kg					
		Main unit: AC adapter (1), p tape fastener roll (1)	ower cable (1), bath guide (1), rear cover (1), single-sided tap	e fastener roll (1), double-sided		
Included accesso	ones	Glass set: Cooling condenser (type A/B)(1), rotary joint (1), evaporation flask (1), receiving flask (1), ball joint clamp (1), flask clip (1), vacuum seal (1), condenser insulation kit (1), condenser bracket (1) (for type B), hex wrench (1) (for type B)					

*1 Performance data above based on 23 ±5 °C room temperature, 65%RH ±20% humidity, and no process load.

*2 Dimensions excludes protrusions.

*3 Applicable rotation speed range and sample volume depend on the capacity of evaporation flask.

Operational Accessories

Glassware Set

Product code	Product code Set					
RG202A	Set A (use with chille Traditional glass set tilted diagonally	,				
RG202B	Set B (use with chille Standard glass set w set vertically, suitable	here condenser is				
Conde	enser A	Condenser B				

Bath Specifications

Product name	Water bath		Oil bath		
Model	BM302-A	BM312-A	BO302-A	BO312-A	
Temp. control range *1	RT +10°C~90°C		RT +10°C~180°C		
Temp. control accuracy *1	±1.0°C		±1.5°C (water), ±2°C (oil)		
Safety features	Automatic overheat prevention, independent overheat prevention (fixed temp.), temperature upper limit error, overcurrent protection fuse				
Other features Calibration offset, overshoot alert, auto re service outlet (for AC100~115V)			ert, auto resume (s	selectable), 2A	
Bath capacity	~5L				
Internal dimensions	Φ240 x H119 mm				
External dimensions *2	Φ262(max depth 286mm) x H240 mm				
Weight	~4.5kg				
Power supply (fuse capacity)	AC100~115V 10~12A (15A)	AC200~230V 5~6A (10A)	AC100~115V 10~12A (15A)	AC200~230V 5~6A (10A)	

*1 Performance data above based on 220 VAC ±5% supplied power, 23 ±5 °C ambient, 65%RH. ±20% humidity, and no process load. Temp accuracy measured in JTM K05. *2 Dimensions excludes protrusions.



Operational Accessories

Vacuum Controller



Equipped with a color LCD that allows you to judge the operating status by color.

- Suppresses bumping with three operation modes that can be selected according to application.
- Evaporator body and vacuum pump work together with one button
- Wireless connection with specified vacuum pumps for control

Model	VR102S		
Setting range of the degree of vacuum*1	0 to 1013 hPa		
Measurement range of the degree of vacuum	0 to 1100 hPa		
Display	Color LCD (2.3")		
Display items	Measurement/setting of vacuum degree, operation time, status		
Operation mode	Manual (constant operation), Gradient (gradient decompression, constant operation) Auto (gradient decompression, target pressure automatic setting)		
Hold function	Maintains current vacuum degree in the middle of decompression (controlled)		
Pressure unit	mmHg/Torr/hPa/kPa/mbar		
Automatic functions at end of operation	Auto leak, Auto cleaning		
Vacuum control system*2	Based on ON/OFF of vacuum control solenoid valve, or pump rotation speed.		
Safety functions	Communication error, pressure sensor error, memory error, leak error, high pressure error, auto leak at error occurrence		
External dimension*3	86AW x 113D x 83H		
Power supply	24V DC *1 (100-240VAC 1A or less)		
Weight	0.5 kg		
Included accessories	Vacuum line branch joint (O.D. φ4 x φ2 x 700 mm with PTFE tube)		

* A separate optional vacuum control solenoid valve is required. With the combination of vacuum pump N820G and vacuum pump control unit G, vacuum control can be performed by pump rotation speed without vacuum control solenoid valve, However, operation mode is limited.

*² An optional connection cable (multi line of power supply/operation signal) for connecting RE unit is required. When using this unit without connecting to RE units, use optional AC adapter/power cable separately. *3 Dimensions excludes protrusions.

Optional Accessories

Glassware



Evaporating flask Size: 24/40 Product code Ca

255713

255712 255711

255710

255709

255708

Size: 29/38

Receiving flask Size 35/20

Capacity	Product code	Capacity
2L	255719	2L
1L*	255718	1L*
500ml	255717	500ml
300ml	255716	300ml
200ml	255715	200ml
100ml	255714	100ml

*Standard

Product code	Capacity
255706	2L
255705	1L
255704	500ml
255703	300ml
255702	200ml
255701	100ml



Rotary joint RG202A					
Product code	Size	Туре			
255720	29/38 L284mm	Standard			
255722	24/40 L286mm	Standard*			
255724	29/38 L284mm	Transparent			
255726	24/40 L286mm	Transparent			
RG202B					
Product code	Size	Туре			
255721	29/38 L208mm	Standard			
255723	24/40 L210mm	Standard*			
255725	29/38 L208mm	Transparent			
255727	24/40 L210mm	Transparent			
	RE202	/212 Series			



Rotary joint different diameter

Product code	Description
255732	$24/40 \rightarrow 24/40 \text{ L105mm}$
255733	$24/40 \rightarrow 19/38 \text{ L103mm}$
255734	$24/40 \rightarrow 15/25 \text{ L90mm}$
255728	$29/38 \rightarrow 29/38 \text{ L106mm}$
255729	29/38 \rightarrow 24/40 L108mm
255731	$29/38 \rightarrow 15/25 \text{ L105mm}$
255730	29/38 → 19/38 L105mm



PTFE vacuum seal					
Product code	Description				
255741	Recommended for ketone and other solvents				



FKM vacuum seal Product code Description RE50040090 YSA Standard

Optional Accessories

Water Circulator (Chiller) CF302LA/CF312LA



Control panel

Model	CF302L-A / CF312L-A	CF802-A		
Operating temp. range	-20°C~30°C			
Temp. control accuracy	±1.0°C (≥ 0°C) ±1.5°C (< 0°C)	±1.0°C		
Cooling capacity	~420W at liquid temp10°C ~280W at liquid temp-10°C	~1320W (at 10°C) ~700W (at -10°C)		
Temp. control	Refrigerator On/Off			
Refrigerator, coolant	Air cooling 450W, R404A	Air cooling, 700W, R410A		
External dimension WxDxH (including protrusions)	205 × 396 × 535 mm (225 × 434 × 564 mm)	340 × 370 × 838 mm (340 × 408 × 920 mm)		
Water capacity	~3.9L (Liquid volume 3.5L)	~15.5L (liquid volume 14L)		
Power source 60Hz	115V 7.5A / 220V 4A	Single phase 115V 15A		
Weight	~31kg	~44kg		



CF802A





Vacuum Pump



Model	N820G
Ultimate vacuum (mbar abs.)	
Minimum speed:	
Gas ballast closed	≤ 6
Gas ballast open	≤ 17
Maximum speed:	
Gas ballast closed	≤ 8
Gas ballast open	≤ 15
Flow rate at atm. pressure (I/min)	
Min. speed	9
Max. speed	20
Permissible ambient temperature	+10°C to + 40°C
Diaphragm material	PTFE-coated
Voltage (V)	100-240
Frequency (Hz)	50/60
Max. operating current (A)	0.66
Dimensions LxWxH (mm)	259 x 163 x 220
Weight (kg)	8.8

Filter mounting plate

Vacuum Pump Control Unit G

255783 (OVR26)

- Regulates motor speed of vacuum pump to control the degree of vacuum
- Wireless interconnection with VR102S

Included accessories:

- N820G mounting bracket 1 pc.
- Hexagon wrench 1 pc.
- M4 hexagon socket head cap screw 3 pcs.

Compatible pump:

Yamato Scientific N820G

Exhaust Trap Kit

255771 (ORT12)



- Used as a solvent recovery unit by installing onto ORT10 stand.
- Comes with an exhaust trap, 500 ml flask, tray, connection hose on the OUT side, and a set of attachment brackets

Vacuum Control Solenoid Valve



255762 (OVR10)

- Opens / closes to control the degree of vacuum
- Works with VR102S Vacuum Controller by wire
- Can be installed onto RT302 solvent recovery unit

Included accessories:

Circulation hose connection

- PTFE tube $\phi4 \times \phi2 \times 700$ mm 1 pc
- Solenoid valve cable 1 pc



255770 (ORT10)

- Stand for VR102S and vacuum pump
- Power supply to this unit is either RE connection cable or AC adapter
- Includes waste liquid trap bottle (250 ml) 255772 (ORT14)

Solvent Recovery Unit

RT302 (255378)


- Allows efficient solvent recovery by cooling water circulation Connected to the exhaust side of a diaphragm
 - vacuum pump
 - Combination of stand (255770) and exhaust trap kit (255771)






•

Set Guide



SET WITH VACUUM CONTROLLER MAIN UNIT + VACUUM CONTROLLER + BATH + GLASSWARE SET




Basic Set Selection Chart

Model	Glassware		Water Bath		Oil Bath		Vacuum
wouer	A	В	BM302	BM312	BO302	BO312	Controller
RE-202-AWA	•		•				
RE-202-BWA		•	•				
RE-202-AOA	•				•		
RE-202-BOA		•			•		
RE-212-AW2A	•			•			
RE-212-BW2A		•		•			
RE-212-AO2A	•					•	
RE-212-BO2A		•				•	
REV-202M-AWA	•		•				•
REV-202M-BWA		٠	•				•
REV-202M-AOA	•				•		•
REV-202M-BOA		٠			•		•
REV-212M-AW2A	•			•			•
REV-212M-BW2A		•		•			•
REV-212M-AO2A	•					•	•
REV-212M-BO2A		•				•	•

SET VARIATION


SET 1 BASIC SET Includes a space-saving glassware set B (vertical condenser) and water bath. Also available in glassware set A (diagonal





SET 2 STANDARD SET

Includes a space-saving glassware set B (vertical condenser), water circulator (chiller), water bath and vacuum pump.



SAMPLE VARIATION

SET 3 COMPLETE SET 1

Similar to Set 2 (glassware set B, water circulator, water bath and vacuum pump) but with the addition of a <u>vacuum controller</u>. User friendly and space saving.



Product code	Description
REV202M	Rotary evaporator with vacuum controller
RG202B	Glassware set B (vertical condenser)
BM302	Water bath
CF302LA	Water circulator (chiller)
221581	Insulation hose 2 pcs.
N820G	Diaphragm vacuum pump
255783	Vacuum pump control unit G (OVR26)
255297	Vacuum hose

SET 4 COMPLETE SET 2

Similar to Set 3, this comes with a rotavap with vacuum controller, glassware B, chiller, water bath and vacuum pump. In addition, it considers environment protection and odor measures by improving collection efficiency through a <u>secondary trap.</u>



Product code	Description
REV202M	Rotary evaporator with vacuum controller
RG202B	Glassware set B (vertical condenser)
BM302	Water bath
255770	Stand (ORT10)
255771	Exhaust trap kit (ORT12)
CF302LA	Water circulator (chiller)
221581	Insulation hose 2 pcs.
281330	Insulation hose (OA094) 1 pc.
N820G	Diaphragm vacuum pump
255783	Vacuum pump control unit G (OVR26)
255297	Vacuum hose

SET 5 "USE YOUR OWN VACUUM PUMP" SET

A complete set with a rotavap with vacuum controller, spacesaving glassware set B (vertical condenser), water circulator (chiller) and water bath, combined with your <u>existing vacuum</u> <u>pump</u>. Vacuum pump is manually and continuously operated.



Product code	Description
REV202M	Rotary evaporator with vacuum controller
RG202B	Glassware set B (vertical condenser)
BM302	Water bath
CF302LA	Water circulator (chiller)
221581	Insulation hose 2 pcs.

SET 6 DOUBLE PERFORMANCE SET

A combination that connects two complete sets of rotary evaporator (glassware set B, water bath and vacuum pump) with one water circulator (CF802) that can be installed under the laboratory table or fume hood, and a secondary trap.



Product code	Qty	Description
REV202M	2	Rotary evaporator with vacuum controller
RG202B	2	Glassware set B (vertical condenser)
BM302	2	Water bath
CF802	1	Water circulator (chiller)
281478	1	Secondary trap (OCF84)
221581	2 sets	Insulation hose 2 pcs.
N820G	2	Diaphragm vacuum pump
255783	2	Vacuum pump control unit G (OVR26)
255297	2	Vacuum hose