

Purchase a Thermo Scientific Orion 3-Star pH Meter from Daigger® and Receive a FREE Replacement Electrode

Purchase a 3-Star Meter



Receive your FREE Electrode

— a \$238.00 value!



Get the best selling meter, the 3-Star pH meter (FX4431A), and receive a free replacement electrode (FX4431C)

- IP-54 Splashproof and dustproof benchtop meter
- Backlit display for easy viewing
- SMART STABILITY™ and SMART AVERAGING™ functions for best accuracy and response time
- Store up to 10 password protected methods
- 1000 GLP data log points
- Autosampler capabilities for high through-put sample applications

To receive your **FREE** Orion Epoxy pH/ATC Triode® Electrode (FX4431C) – a \$238.00 value; you **MUST** mention code SFQ2PH1 at the time you place your Orion 3-Star meter order.

FREE replacement electrode will be shipped with order.

Specifications	
Mode	pH, mV, ORP, T
pH Range	0 to 14.00 pH
pH Accuracy	± 0.002 pH
mV Range	± 1999.9
mV Accuracy	± 0.05%
Temp Range	- 5.0 to 105 ° C
Temp Accuracy	± 0.1 ° C
EPA Compliant	CE, UL, CSA
Power Requirement	Universal power adapter
Warranty	3 Years
Ordering Information	
Cat. No.	FX4431A
Description	3-Star pH meter
List Price	\$ 707
Sale Price	\$ 609
Promo Code	SFQ2PH1
Valid Date:	Offer good till June 30,2008

Call Daigger® at 800-621-7193 to take advantage of this offer.

DAIGGER®

P-DAIGGER-E 03/08 RevA

Thermo
SCIENTIFIC