

VAN GUARD

1200 SERIES Clinical Microscopes



Features

- Precision Optics
- 20W Variable Quartz Halogen Illumination
- Brightfield & Phase Contrast
- Integrated Koehler Illumination
- 5-Year Warranty

VAN GUARD[®]
INFINITE • SOLUTIONS



1200-Series Clinical Microscopes


VanGuard 1200-Series Clinical Microscopes provide the necessary specifications for diagnostic & educational use in medical/veterinary clinics, universities, and industrial laboratories.

Wide ranges of models are available in Brightfield & Phase Contrast configurations. Choose from standard achromatic or plan achromatic objectives; monocular, binocular, or trinocular viewing heads.

Various camera packages are available as options; including digital, video, and 35mm systems.

All models feature 10X WideField eyepieces, 20W variable quartz halogen, field diaphragms for Koehler illumination, and are backed by a 5-Year warranty.

Make your lab work better...make it VanGuard.



Technical Specifications

	1210CM	1220CM	1230CM	1211CM	1221CM	1231CM	1212CM	1222CM	1232CM	1213CM	1223CM	1233CM
Optical System	Fixed 160mm			Fixed 160mm			Fixed 160mm			Fixed 160mm		
Type	Monocular	Binocular	Trinocular	Monocular	Binocular	Trinocular	Monocular	Binocular	Trinocular	Monocular	Binocular	Trinocular
Tube Inclination	30°			30°			30°			30°		
Rotatable	360°			360°			360°			360°		
Interpupillary	N/A	55-75mm		N/A	55-75mm		N/A	55-75mm		N/A	55-75mm	
Dioptric	N/A	-5 to +5		N/A	-5 to +5		N/A	-5 to +5		N/A	-5 to +5	
Camera Port	N/A	N/A	20/80 Split	N/A	N/A	20/80 Split	N/A	N/A	20/80 Split	N/A	N/A	20/80 Split
Eyepieces	WideField			WideField			WideField			WideField		
Magnification	10X			10X			10X			10X		
Field Diameter	18mm			18mm			18mm			18mm		
High-Eyepoint	Yes			Yes			Yes			Yes		
Rubber Eyecups	Yes			Yes			Yes			Yes		
Nosepiece	4			4			4			4		
Pitch	Reverse			Reverse			Reverse			Reverse		
Illumination Technique(s)	Brightfield			Brightfield			Brightfield / Phase Contrast			Brightfield / Phase Contrast		
Objectives	Achromatic			PlanAchromatic			Achromatic			PlanAchromatic		
Magnification (N.A.)	4X (0.10), 10X (0.25), 40X (0.60), 100X (1.25)			4X (0.10), 10X (0.25), 40X (0.60), 100X (1.25)			10X (0.25), 20X (0.40), 40X (0.60), 100X (1.25)			10X (0.25), 20X (0.40), 40X (0.60), 100X (1.25)		
DIN Standard	Yes			Yes			Yes			Yes		
Stage	Mechanical			Mechanical			Mechanical			Mechanical		
Shape	Rectangular			Rectangular			Rectangular			Rectangular		
Dimensions	160 x 140mm			160 x 140mm			160 x 140mm			160 x 140mm		
X-Y Control	Drop-Down, Right-Hand, Rack & Pinion			Drop-Down, Right-Hand, Rack & Pinion			Drop-Down, Right-Hand, Rack & Pinion			Drop-Down, Right-Hand, Rack & Pinion		
Cross Travel Range	75 x 50mm			75 x 50mm			75 x 50mm			75 x 50mm		
Slide Holder	Removable Spring Clip			Removable Spring Clip			Removable Spring Clip			Removable Spring Clip		
Finish	Chemical Resistant			Chemical Resistant			Chemical Resistant			Chemical Resistant		
Focusing	Coaxial			Coaxial			Coaxial			Coaxial		
Controls	Coarse / Fine			Coarse / Fine			Coarse / Fine			Coarse / Fine		
Tension Adjustment	Yes			Yes			Yes			Yes		
Safety Auto-Stop	Yes			Yes			Yes			Yes		
Focusing Range	40mm			40mm			40mm			40mm		
Precision	0.002mm			0.002mm			0.002mm			0.002mm		
Condenser	Abbe 1.25 N.A. w/ Iris Diaphragm			Abbe 1.25 N.A. w/ Iris Diaphragm			Abbe 1.25 N.A. w/ Iris Diaphragm			Abbe 1.25 N.A. w/ Iris Diaphragm		
Focusing Control	Rack & Pinion			Rack & Pinion			Rack & Pinion			Rack & Pinion		
Zernike Phase Annuli	N/A			N/A			10X, 20X, 40X, & 100X			10X, 20X, 40X, & 100X		
Illumination	20W Variable Quartz Halogen			20W Variable Quartz Halogen			20W Variable Quartz Halogen			20W Variable Quartz Halogen		
Type	Field Diaphragm			Field Diaphragm			Field Diaphragm			Field Diaphragm		
Koehler	On Collector			On Collector			On Collector			On Collector		
Filter Holder	On Collector			On Collector			On Collector			On Collector		
Base	210 x 180mm; Anti-Skid Rubber Feet			210 x 180mm; Anti-Skid Rubber Feet			210 x 180mm; Anti-Skid Rubber Feet			210 x 180mm; Anti-Skid Rubber Feet		
Body	Ergonomic; Stain-Resistant Enamel Finish			Ergonomic; Stain-Resistant Enamel Finish			Ergonomic; Stain-Resistant Enamel Finish			Ergonomic; Stain-Resistant Enamel Finish		
Dimensions	210 (L) x 180 (W) x 380 (H) mm			210 (L) x 180 (W) x 380 (H) mm			210 (L) x 180 (W) x 380 (H) mm			210 (L) x 180 (W) x 380 (H) mm		
Weight	7.0kg			7.0kg			7.0kg			7.0kg		
Voltage	110V (220V Available On Request)			110V (220V Available On Request)			110V (220V Available On Request)			110V (220V Available On Request)		
Warranty	5-Year Optical/Mechanical; 1-Year Electrical			5-Year Optical/Mechanical; 1-Year Electrical			5-Year Optical/Mechanical; 1-Year Electrical			5-Year Optical/Mechanical; 1-Year Electrical		
Supplied With	Clear Diffusing Filter, Blue Filter, Spare Lamps & Fuses, Dust Cover, & Operation Manual			Clear Diffusing Filter, Blue Filter, Spare Lamps & Fuses, Dust Cover, & Operation Manual			Clear Diffusing Filter, Blue Filter, Spare Lamps & Fuses, Telescoping Eyepiece, Spare Lamps & Fuses, Dust Cover, & Operation Manual			Clear Diffusing Filter, Blue Filter, Spare Lamps & Fuses, Telescoping Eyepiece, Spare Lamps & Fuses, Dust Cover, & Operation Manual		



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
 www.veege.com • e-mail: sales@veege.com

WARNING: This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov