



Volumetric Flasks

Class A, Φ Stopper, To Contain

SIBATA 2306A-Series Volumetric Flasks meet ASTM E-288 Class A Specifications & USP Standards for Volumetric Glassware.

These flasks feature machine-blown bodies and heavy-beaded, heavy-tubing necks, tooled for Φ stoppers. The graduation line is sharp and permanent; large, white, block-letters allow the labeling to be easily read. Each flask also has a large marking square.


All sizes are supplied with a glass Φ stopper.

Specifications

Cat. No.	Capacity (mL)	Tolerance (\pm mL)	Φ Stopper No.	Approx. O.D. x Height (mm)
2306A-5	5	0.02	9	21 x 72
2306A-10	10	0.02	9	26 x 96
2306A-25	25	0.03	9	37 x 123
2306A-50	50	0.05	13	49 x 152
2306A-100	100	0.08	13	58 x 181
2306A-200	200	0.10	16	74 x 220
2306A-250	250	0.12	16	78 x 239
2306A-500	500	0.20	19	96 x 279
2306A-1000	1000	0.30	22	122 x 327
2306A-2000	2000	0.50	27	157 x 396

Packaging

Cat. No.	Pack Qty.	Case Qty.
2306A-5	6	72
2306A-10	6	72
2306A-25	6	72
2306A-50	4	48
2306A-100	4	48
2306A-200	4	24
2306A-250	4	24
2306A-500	2	12
2306A-1000	2	8
2306A-2000	2	8



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
 www.veege.com • e-mail: sales@veege.net

Printed in the USA

WARNING: This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov

VGBCR2306A-Series08292018