

VEE GEE®

Burets

Class A, w/ PTFE Stopcocks

VEE GEE 20114-Series Burets meet ASTM E-287 Class A Specifications.

Made from precision, uniform Schott DURAN® true-bore tubing to ensure the stipulated capacity tolerances. These burets feature durable, permanent markings, fine, sharp lines, and large, easy-to-read numbers. Supplied with PTFE stopcocks to eliminate freezing and lubricant contamination and to provide a precision fit within the barrel.

All sizes are supplied with a 2mm stopcock bore plug.

Specifications

Cat. No.	Capacity (mL)	Graduation Interval (mL)	Tolerance (\pm mL)	Stopcock Bore (mm)	Approx. O.D. x Height (mm)
20114-10	10	0.05	0.02	2	9 x 510
20114-25	25	0.10	0.03	2	12 x 550
20114-50	50	0.10	0.05	2	13 x 758
20114-100	100	0.20	0.10	2	17 x 775

Packaging

Cat. No.	Pack Qty.	Case Qty.
20114-10	1	24
20114-25	1	24
20114-50	1	24
20114-100	1	24

⚠ WARNING: This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veegee.com • e-mail: sales@veegee.com