



# Graduated Cylinders

## Class A, Double Metric Scale, Graduated, To Contain

SIBATA 2351A-Series Graduated Cylinders meet ASTM E-1272 Class A, Style I Specifications & USP Standards for Volumetric Glassware.

Feature hexagonal bases for added stability. All cylinders, with the exception of the 2000mL size, come with protective, polyethylene, hexagonal bumper guards.


All markings are in durable, white enamel.

### Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Tolerance ( $\pm$ mL)	Approx. O.D. x Height (mm)
2351A-5	5	0.5 to 5.0	0.1	0.05	7 x 126
2351A-10	10	0.7 to 10.0	0.1	0.05	8 x 182
2351A-25	25	1.4 to 25.0	0.2	0.15	17 x 210
2351A-50	50	3.5 to 50.0	0.5	0.25	24 x 218
2351A-100	100	5.0 to 100.0	1.0	0.40	30 x 236
2351A-250	250	10.0 to 250.0	2.0	1.00	41 x 301
2351A-500	500	25.0 to 500.0	5.0	2.00	53 x 362
2351A-1000	1000	50.0 to 1000.0	10.0	3.00	70 x 422
2351A-2000	2000	100.0 to 2000.0	20.0	6.00	90 x 478

### Packaging

Cat. No.	Pack Qty.	Case Qty.
2351A-5	4	48
2351A-10	4	48
2351A-25	4	48
2351A-50	4	48
2351A-100	4	48
2351A-250	4	24
2351A-500	2	12
2351A-1000	2	8
2351A-2000	1	6



## VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034  
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826  
 www.veege.com • e-mail: sales@veege.com

Printed in the USA

**WARNING:** This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

VGBRC2351A-Series08292018