



Graduated Cylinders

Class B, Double Metric Scale, Graduated, To Contain

SIBATA 2351-Series Graduated Cylinders meet ASTM E-1272 Class B, Style I Specifications.

Feature hexagonal bases for added stability. All cylinders, with the exception of the 2000mL size, come with protective, polyethylene, hexagonal bumper guards.

All markings are in durable, white enamel.


Specifications

Cat. No.	Capacity (mL)	Graduation Range (mL)	Graduation Interval (mL)	Tolerance (\pm mL)	Approx. O.D. x Height (mm)
2351-10	10	0.7 to 10.0	0.1	0.10	8 x 182
2351-25	25	1.4 to 25.0	0.2	0.30	17 x 210
2351-50	50	3.5 to 50.0	0.5	0.40	24 x 218
2351-100	100	5.0 to 100.0	1.0	0.60	30 x 236
2351-250	250	10.0 to 250.0	2.0	1.40	41 x 301
2351-500	500	25.0 to 500.0	5.0	2.60	53 x 362
2351-1000	1000	50.0 to 1000.0	10.0	5.00	70 x 422
2351-2000	2000	100.0 to 2000.0	20.0	10.00	90 x 478

Packaging

Cat. No.	Pack Qty.	Case Qty.
2351-10	4	48
2351-25	4	48
2351-50	4	48
2351-100	4	48
2351-250	4	24
2351-500	2	12
2351-1000	2	8
2351-2000	1	6

⚠ WARNING: This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veegee.com • e-mail: sales@veegee.com