



Media Bottles

VEE GEE Media Bottles are commonly used to contain chemicals in liquid or powder form in laboratories. The bottles are made from borosilicate glass in compliance with ASTM 438 Type 1 and are available in clear or amber (low actinic) for light sensitive chemical compounds.

Five sizes available in both clear and amber, VEE GEE media bottles have threaded on Blue Plastic caps and pouring rings that are autoclavable.

- Clear and Amber
- Made from Borosilicate Glass
- Threaded Blue Plastic Caps
- Autoclavable Caps & Pouring Ring
- Clear White Marking

Media Bottles, Clear

Cat. No.	Capacity (mL)	Pack Qty.
20274-100	100	10
20274-250	250	10
20274-500	500	10
20274-1000	1000	10
20274-2000	2000	4

Media Bottles, Amber

Cat. No.	Capacity (mL)	Pack Qty.
20275-100	100	10
20275-250	250	10
20275-500	500	10
20275-1000	1000	10
20275-2000	2000	4



WARNING: This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veegee.com • email: sales@veegee.com

VGBRC082918/20270-Series
Printed in the USA