



Serological Pipets

Class A, Reusable Glass, Color-Coded, w/ Colored Markings, To Deliver

SIBATA 2011A-Series Serological Pipets meet ASTM E-1044 Style 1 Specifications.

These pipets are graduated to the tip (tempered) and are also calibrated to deliver. Feature permanent amber markings and color-codings which are enameled onto the glass.

The top end of sizes 5mL to 25mL is constricted.

Specifications

Cat. No.	Capacity (mL)	Tolerance (±mL)	Graduation Interval (mL)	Approx. O.D. x Height (mm)	Color Code
2011A-01-c	0.1	0.0025	0.01	6 x 300	White
2011A-02-c	0.2	0.004	0.01	6 x 300	Black
2011A-1-d	1.0	0.01	0.1	7 x 360	Red
2011A-1-c	1.0	0.01	0.01	7 x 360	Yellow
2011A-2-d	2.0	0.01	0.1	7 x 360	Green
2011A-5	5.0	0.02	0.1	8 x 365	Blue
2011A-10	10.0	0.03	0.1	10 x 370	Orange
2011A-25	25.0	0.05	0.1	13 x 455	White

Packaging

Cat. No.	Pack Qty.	Case Qty.
2011A-01-c	12	72
2011A-02-c	12	72
2011A-1-d	12	72
2011A-1-c	12	72
2011A-2-d	12	72
2011A-5	12	72
2011A-10	12	72
2011A-25	12	72

WARNING: This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veege.com • e-mail: sales@veege.net