



Volumetric Pipets

Class A, Reusable Glass, Color-Coded, w/ Colored Markings, To Deliver

SIBATA 2040A-Series Volumetric Pipets meet ASTM E-969 Class A Specifications & USP Standards for Volumetric Glassware.

These pipets have sharp and permanent graduation lines and are calibrated to deliver. Feature permanent amber markings and color-codings which are enameled onto the glass.



Specifications

Cat. No.	Capacity (mL)	Tolerance (\pm mL)	Approx. Length (mm)	Color Code	Pack Qty.	Case Qty.
2040A-05	0.5	0.005	300	2-Black	12	72
2040A-1	1.0	0.006	345	Blue	12	72
2040A-2	2.0	0.006	355	Orange	12	72
2040A-3	3.0	0.01	360	Black	12	72
2040A-4	4.0	0.01	365	2-Red	12	72
2040A-5	5.0	0.01	375	White	12	72
2040A-6	6.0	0.01	390	2-Orange	12	72
2040A-7	7.0	0.01	410	2-Green	12	72
2040A-8	8.0	0.02	430	Blue	12	72
2040A-9	9.0	0.02	436	Black	12	72
2040A-10	10.0	0.02	440	Red	12	72
2040A-13	13.0	0.03	460	White	6	36
2040A-15	15.0	0.03	460	Green	6	36
2040A-20	20.0	0.03	475	Yellow	6	36
2040A-25	25.0	0.03	480	Blue	6	36
2040A-50	50.0	0.05	540	Red	6	36
2040A-100	100.0	0.08	600	Yellow	6	36

WARNING: This product can expose you to chemicals including lead, which is known to the state of California to cause cancer and birth defects or other reproductive harm. For more information go to: www.P65Warnings.ca.gov



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
(800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
www.veegee.com • e-mail: sales@veegee.net