



Wide-Tip Volumetric Pipets

Reusable Glass, To Deliver

SIBATA 2040-Series Wide-Tip Volumetric Pipets meet ASTM E-969 Class B Specifications for Tolerance.

The wide-tip opening facilitates rapid fill and drain of viscous liquids, slurries, and suspensions; an excellent product for water treatment facilities, winemakers, and essential oil testing.

These pipets have sharp and permanent graduation lines and are calibrated to deliver. Feature permanent amber markings.


Sizes 1, 2, & 5mL do not have drawn tips, while sizes 10mL through 100mL feature standard drawn tips.

Specifications

Cat. No.	Capacity (mL)	Tolerance (\pm mL)	Approx. Length (mm)	Approx. Tip I.D. (mm)
2040-1-1.9MM	1.0	0.01	320	1.9
2040-2-2.4MM	2.0	0.01	330	2.4
2040-5-3MM	5.0	0.02	346	3.0
2040-10-3MM	10.0	0.02	440	3.0
2040-20-3MM	20.0	0.03	478	3.0
2040-25-3MM	25.0	0.03	475	3.0
2040-50-3MM	50.0	0.05	538	3.0
2040-100-3MM	100.0	0.10	605	3.0

Packaging

Cat. No.	Pack Qty.	Case Qty.
2040-1-1.9MM	12	72
2040-2-2.4MM	12	72
2040-5-3MM	12	72
2040-10-3MM	12	72
2040-20-3MM	6	36
2040-25-3MM	6	36
2040-50-3MM	6	36
2040-100-3MM	6	36



VEE GEE Scientific, Inc.

13600 NE 126th Place, Suite A • Kirkland, WA 98034
 (800) 423-8842 • (425) 823-4518 • Fax: (425) 820-9826
 www.veegee.com • e-mail: sales@veegee.net